

Products to be used for  
"Human Health"
























**AILEE**

**SUTURES**

# Contents

1.	NEEDLE TYPE	2
2.	EXPLANATION SUTURE CODE	3
3.	PRODUCT IMAGE	4
4.	SURGIFIT (WHITE, VIOLET)	5
5.	MONOFIT-S (UNDYED, VIOLET)	17
6.	MONOFIT-L (UNDYED, VIOLET)	29
7.	CHROMIC	41
8.	SILK (WHITE, BLACK)	49
9.	NYLON (WHITE, BLUE, BLACK)	65
10.	POLYESTER (WHITE, GREEN)	81
11.	POLYPROPYLENE	89
12.	SUTURE NEEDLES	101
13.	SURGICAL BLADES	113


# NEEDLE TYPE

Type	Shape	Design	Cutting Pattern
TAPER POINT			
BLUNT POINT			
TAPERCUT			
REVERSE CUTTING			
SPATULA			
DIAMOND			
TROCAR POINT			

Straight	Half Curved	1/4 Circle	3/8 Circle
----------	-------------	------------	------------



1/2 Circle	5/8 Circle	Fish hook	Compound Curved
------------	------------	-----------	-----------------



REGULAR EYE	FRENCH EYE	FRENCH SPRING EYE	EYELESS
-------------	------------	-------------------	---------



A



B



C



D

# EXPLANATION SUTURE CODE



**S**

## Denotes the Kind of suture material

- C** = CHROMIC catgut
- A** = Synthetic Absorbable Suture braided
- M** = Synthetic Absorbable Suture Monofilament
- S** = SILK braided and coated with silicone
- N** = NYLON monofilament
- PE** = POLYESTER braided and coated with silicone
- PP** = POLYPROPYLENE monofilament

**W**

## Denotes the color of suture

- W** = White
- B** = Blue
- G** = Green
- V** = Violet
- K** = Black

**6**

## Denotes the suture size in USP (EP)

- |                           |                        |
|---------------------------|------------------------|
| <b>11</b> = 11-0 (EP 0.1) | <b>3</b> = 3-0 (EP 2 ) |
| <b>10</b> = 10-0 (EP 0.2) | <b>2</b> = 2-0 (EP 3 ) |
| <b>9</b> = 9-0 (EP 0.3)   | <b>1</b> = 0 (EP 3.5)  |
| <b>8</b> = 8-0 (EP 0.4)   | <b>01</b> = 1 (EP 4 )  |
| <b>7</b> = 7-0 (EP 0.5)   | <b>02</b> = 2 (EP 5 )  |
| <b>6</b> = 6-0 (EP 0.7)   | <b>03</b> = 3 (EP 6 )  |
| <b>5</b> = 5-0 (EP 1 )    | <b>04</b> = 4 (EP 7 )  |
| <b>4</b> = 4-0 (EP 1.5)   | <b>05</b> = 5 (EP 7 )  |

**10**

## Denotes the serial number

(or suture length in case non-needed sutures)

**DN**

## Denotes speciality of needle & suture

- D** =Diamond shape needle
- DN**=Double needle
- H** =Heavy needle
- L** =Loop needle
- M** =Multi-Strand suture
- P** =Plastic Surgery needle
- R** =Round shape needle
- RM**=Removable needle
- T** =Tapercut shape needle



# SURGITIT

(Braided Polyglycolic acid sutures)

## Surgift (써지피트 브랜드 흡수성봉합사)

PGA (Polyglycolide Suture)

### DESCRIPTION 설명

100% Polyglycolide 합성사로 구성된 흡수성봉합사로 가수분해 되어 지며, 60일에서 90일 안에 완전히 흡수된다. 복선사로 되어있어 핸들링이 용이하여 매듭의 안전성이 뛰어나다.

Poly(caprolactone)와 Calcium Stearate로 코팅되어 조직을 부드럽게 통과하고 쉽게 매듭이 내려가도록 재질의 상호작용을 도와주고, 높은 초기 인장강도로 상처가 치유되는 기간 동안 강한 힘으로 지탱해주며 그 후에는 빠르게 흡수된다.

100% synthetic absorbable braided Surgifit Suture, which is made from P.G.A(Polyglycolic Acid), can be hydrolyzed and fully absorbable within 60 days ~90 days and it gives easy passage-through and superior knot stability. It is coated with Poly(caprolactone) and Calcium Stearate, so it ensures easy passage -through, and its very high initial tensile strength preserves the stitched part strongly during healing period, and later it shows fast absorption.

### COLOUR 색상

Violet, Beige(Undyed)

### STERILIZATION 멸균

Ethylene Oxide

### SIZE 사이즈

USP 8/0 ~ 2

EP 0.4 ~ 5

% of original strength remaining

2 Weeks	75%
3 Weeks	50%
4 Weeks	25%

### CERTIFICATION




- Certificate of GMP (KFDA)
- ISO9001
- ISO13485
- CE Mark

















※ Catalogue에 명시된 제품 외에도 주문에 따라 생산 가능 합니다.

# SURGITIT

(Braided Polyglycolic acid sutures)

SURGITIT

Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	8-0 (0.4)	7-0 (0.5)	6-0 (0.7)	5-0 (1)
<b>6m/m</b>								
3/8 	 220 Microns	Violet	30cm	AV8221DN	AV7221DN			
<b>8m/m</b>								
1/4 	 220 Microns	Violet	30cm	AV825P				
1/4 		White	45cm				AW6081DN	AW5081DN
		Violet	45cm				AV6081DN	AV5081DN
3/8 	 220 Microns	Violet	30cm				AV6073	
1/2 		White	45cm				AW6072DN	
1/2 	 220 Microns	Violet	45cm			AV7262DN		
<b>9m/m</b>								
3/8 		White	45cm				AW6101	AW5101
		Violet	45cm				AV6101	







Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	7-0 (0.5)	6-0 (0.7)	5-0 (1)	4-0 (1.5)
<b>9m/m</b>								
3/8			Violet	45cm			AV6101DN	
3/8			White	45cm			AW510DN	
1/2			Violet	45cm			AV611	
<b>11 m/m</b>								
3/8			White	45cm			AW617P	AW517P
			Violet	45cm			AV617P	AV517P
<b>13m/m</b>								
			Violet	45cm				AV421S
3/8				45cm			AV621R	
			Violet	60cm				AV421R
				70cm			AV521R	
3/8			White	45cm			AW621	AW521
			Violet	45cm			AV621	AV521
1/2			White	45cm			AW6220	AW5220
			Violet	45cm				AV5220















# SURGITIT

(Braided Polyglycolic acid sutures)







SURGITIT











Needle		Suture		Product Code					
Length	Description	Color	Length	USP (EP)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)
<b>14m/m</b>									
1/2		White	45cm		AW622	AW522	AW422		
		Violet	45cm		AV622	AV522	AV422		
		Violet	70cm		AV6221	AV5221	AV4221		
<b>15m/m</b>									
3/8		White	45cm					AW304P	
3/8		White	45cm			AW526	AW426		
1/2		Violet	70cm			AV5261			
<b>17m/m</b>									
3/8		Violet	45cm				AV4291		
3/8		White	70cm			AW528	AW428		
		Violet	70cm			AV528	AV428	AV306	
1/2		Violet	70cm			AV529	AV429	AV307	AV207

Needle		Suture		Product Code				
Length	Description	Color	Lenght	USP (EP)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)
<b>18m/m</b>								
3/8		Violet	70cm		AV530	AV430	AV308	
								
3/8		White	45cm		AW534	AW434		
				70cm				AW312
				45cm			AV434	
<b>19m/m</b>								
3/8		White	45cm			AW4341	AW3121	
				70cm			AW4342	AW3122
				45cm				AV3121
1/2		White	70cm			AW435		
				70cm		AV535	AV435	AV313
<b>21 m/m</b>								
3/8		Violet	45cm			AV441		
								
3/8		Violet	70cm				AV317	AV217
								

# SURGITIT








(Braided Polyglycolic acid sutures)







Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	4-0 (1.5)	3-0 (2)	2-0 (3)	1-0 (3.5)
<b>21 m/m</b>								
1/2	⊙	Violet	70cm	AV4391	AV3171	AV2171		
								
1/2	⊕	Violet	70cm	AV402T				
								
1/2	▽	Violet	45cm		AV320			
								
<b>24 m/m</b>								
1/2	⊙	White	70cm		AW322			
		Violet	70cm	AV444	AV322			
<b>25 m/m</b>								
3/8	▽	White	70cm	AW446P	AW324	AW224		
		Violet	70cm		AV324	AV224		
1/2	⊙	White	70cm	AW447				
		Violet	70cm	AV447 AV447P	AV325 AV325P	AV225 AV225P		AV107

Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	3-0 (2)	2-0 (3)	1-0 (3.5)	1 (4)
<b>28 m/m</b>								
1/2								
		Violet	70cm		AV326	AV226	AV108	AV0108
5/8								
		Violet	70cm		AV326E	AV226E	AV126E	
<b>29 m/m</b>								
1/2								
		White	70cm			AW229		
		Violet	70cm		AV329	AV229	AV111	AV0111
<b>31 m/m</b>								
1/2								
		Violet	70cm			AV230	AV112H AV112	AV0112
<b>33 m/m</b>								
1/2								
		Violet	70cm			AV235		
		Violet	90cm				AV1171	

# SURGITIT

(Braided Polyglycolic acid sutures)





Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	3-0 (2)	2-0 (3)	1-0 (3.5)	1 (4)
<b>36 m/m</b>								
1/2		White	90cm			AW236		
		Violet	70cm			AV236P	AV118	
		Violet	90cm			AV236		AV0118
<b>37 m/m</b>								
3/8		Violet	70cm			AV239		
1/2		Violet	70cm	AV340		AV240	AV122	AV0122
1/2		Violet	70cm				AV122T	
<b>41 m/m</b>								
1/2			70cm			AV244	AV126	AV0127
		Violet	75cm				AV127	
		Violet	90cm			AV2271	AV1271	AV01271

Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	2-0 (3)	1-0 (3.5)	1 (4)	2 (5)
<b>41 m/m</b>								
1/2	☉							
		Violet	70cm				AV0127T	
1/2	▽							
		Violet	70cm			AV124	AV0124	
<b>45 m/m</b>								
3/8	▽							
		Violet	70cm			AV131		
1/2	☉							
		Violet	90cm				AV0132	AV0232
<b>46 m/m</b>								
1/2	☉							
		Violet	70cm			AV133	AV0133	AV0233
1/2	☉							
		Violet	70cm				AV0133T	AV0233T

# SURGITIT

(Braided Polyglycolic acid sutures)

SURGITIT

Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	2-0 (3)	1-0 (3.5)	1 (4)	2 (5)
<b>49 m/m</b>								
1/2								
		Violet	70cm				AV0151	AV0251
			90cm				AV01511	
1/2								
		Violet	90cm					AV0252

## Non-needled Suture

Suture		Product Code										
Color	Length	USP (EP)	8-0 (0.4)	7-0 (0.5)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)	1-0 (3.5)	1 (4)	2 (5)
White Violet	45cm x 1		AV 8451	AV 7451	AW 6451 AV 6451	AW 5451 AV 5451	AW 4451 AV 4451	AW 3451 AV 3451	AW 2451 AV 2451	AV 1451	AV 01451	AV 02451
White Violet	100cm x 1		AV 81001	AV 71001	AW 61001 AV 61001	AW 51001 AV 51001	AW 41001 AV 41001	AW 31001 AV 31001	AW 21001 AV 21001	AV 11001	AV 011001	AV 021001
White Violet	120cm x 1		AV 81201	AV 71201	AW 61201 AV 61201	AW 51201 AV 51201	AW 41201 AV 41201	AW 31201 AV 31201	AW 21201 AV 21201	AV 11201	AV 011201	AV 021201
White Violet	135cm x 1					AW 51351 AV 51351	AW 41351 AV 41351	AW 31351 AV 31351	AW 21351 AV 21351	AV 11351	AV 011351	AV 021351
White Violet	150cm x 1		AV 81501	AV 71501	AW 61501 AV 61501	AW 51501 AV 51501	AW 41501 AV 41501	AW 31501 AV 31501	AW 21501 AV 21501	AV 11501	AV 011501	AV 021501
White Violet	180cm x 1		AV 81801	AV 71801	AW 61801 AV 61801	AW 51801 AV 51801	AW 41801 AV 41801	AW 31801 AV 31801	AW 21801 AV 21801	AV 11801	AV 011801	AV 021801





# MONOFIT-S

(Monofilament  
Polyglycolic-co-caprolactone  
sutures)

## Monofit-S (모노피트-S 흡수성봉합사)

PGCL (Polyglycolic-co-caprolactone Suture)

### DESCRIPTION 설명

75% Glycolide와 25% Caprolactone의 중합체로 구성된 흡수성합성사로 가수분해 되어 지며, 80일에서 120일 안에 완전히 흡수된다.

한 가닥으로 구성된 단선사로서 표면이 부드러워 조직을 부드럽게 통과하고 조직 통과 시 낮은 마찰력을 가지고 있다. 또한 우수한 유연성으로 실의 꼬임현상이 적다. 또한 높은 초기 인장강도로 상처가 치유되는 기간 동안 강한 힘으로 지탱해주고 그 후에는 빠르게 흡수된다.

Polymer absorbable suture, which is made from 75% Glycolide and 25% Caprolactone, can be hydrolyzed and fully absorbable within 80 ~ 120 days and it gives easy passage-through with its single strand suture and it has low frictional force during tissue passage-through. And, its superior flexibility gives low suture twist. And, its very high initial tensile strength preserves the stitched part strongly during healing period, and later it shows fast absorption.

### COLOUR 색상

Violet, Beige(Undyed)

### STERILIZATION 멸균

Ethylene Oxide

### SIZE 사이즈

USP 6/0 ~ 1

EP 0.7 ~ 4

% of original strength remaining

1 Weeks	50~55%
2 Weeks	20%
3 Weeks	5%

### CERTIFICATION















- Certificate of GMP (KFDA)
- ISO9001
- ISO13485

※Catalogue에 명시된 제품 외에도 주문에 따라 생산 가능 합니다.

# MONOFIT-S

(Monofilament Polyglycolic-co-caprolactone sutures)

MONOFIT-S

Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)
<b>9m/m</b>								
3/8 		Violet	45cm	MVS6012	MVS510			
1/2 		Violet	45cm	MVS6011	MVS511			
<b>11 m/m</b>								
		Violet	45cm				MVS4131DN	
3/8 		Undyed	45cm	MWS617	MWS517	MWS417		
		Violet	45cm	MVS617	MVS517	MVS417		
3/8 		Violet	45cm		MVS5012			
3/8 		Undyed	45cm	MWS602	MWS501			
		Violet	45cm	MVS602	MVS501			
1/2 		Undyed	45cm	MWS6013				
		Violet	45cm	MVS6013	MVS513	MVS413		

# MONOFIT-S

(Monofilament Polyglycolic-co-caprolactone sutures)






MONOFIT-S

Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)
<b>13m/m</b>								
3/8	⊙	Violet	70cm				MVS402	
	☾							
3/8	▽	Undyed	45cm	MWS621	MWS5021			
		Violet	45cm	MVS621	MVS5021			
3/8	△	Violet	45cm			MVS5022		
3/8	◻	Undyed	45cm	MWS603	MWS502	MWS401		
		Violet	45cm	MVS603	MVS502	MVS401		
1/2	⊙	Violet	45cm	MVS604	MVS504			
			70cm		MVS503			
1/2	▽	Violet	45cm	MVS619	MVS519			

# MONOFIT-S

(Monofilament Polyglycolic-co-caprolactone sutures)









MONOFIT-S

Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)
<b>14m/m</b>								
1/2		Violet	45cm	MVS622	MVS522	MVS422		
								
<b>15m/m</b>								
3/8		Violet	70cm				MVS426R	
3/8		Undyed	70cm				MWS426	
				45cm	MVS626	MVS526	MVS4031	
		Violet	70cm				MVS426	
3/8		Violet	45cm				MVS4032	

# MONOFIT-S

(Monofilament Polyglycolic-co-caprolactone sutures)









MONOFIT-S

Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)
<b>17 m/m</b>								
3/8								
		Violet	70cm				MVS405	
3/8								
		Undyed	70cm				MWS428	MWS3011
			70cm				MVS428	MVS3011
		Violet	75cm			MVS528		
3/8								
		Violet	45cm			MVS505		
1/2								
		Undyed	45cm			MWS529	MWS429	
			45cm			MVS529	MVS429	
		Violet	70cm				MVS406	MVS301

# MONOFIT-S

(Monofilament Polyglycolic-co-caprolactone sutures)







MONOFIT-S

Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)
<b>18m/m</b>								
3/8 		Undyed	45cm				MWS434 MWS404	MWS312
		Violet	45cm				MVS434 MVS404	MVS312 MVS303
3/8 		Violet	70cm			MVS506		
1/2 		Undyed	45cm				MWS430	
		Violet	45cm	MVS630	MVS530	MVS430		
<b>19m/m</b>								
1/2 		Violet	45cm				MVS435	MVS313

# MONOFIT-S

(Monofilament Polyglycolic-co-caprolactone sutures)

MONOFIT-S







Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)
<b>21 m/m</b>								
3/8 		Undyed	45cm			MWS439		
		Violet	45cm			MVS439	MVS317	
1/2 		Undyed	70cm			MWS442	MWS320	
		Violet	70cm			MVS442	MVS320	MVS220
<b>23 m/m</b>								
1/2 		Violet	70cm			MVS443	MVS321	MVS221



# MONOFIT-S

(Monofilament Polyglycolic-co-caprolactone sutures)









MONOFIT-S

Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	4-0 (1.5)	3-0 (2)	2-0 (3)	1-0 (3.5)
<b>25 m/m</b>								
3/8 		Violet	45cm	MVS446				
			70cm		MVS306	MVS204		
1/2 		Undyed	70cm		MWS325	MWS205		
			70cm	Violet	MVS447	MVS325	MVS205	
<b>28 m/m</b>								
1/2 		Violet	70cm		MVS326	MVS226	MVS101	

# MONOFIT-S

(Monofilament Polyglycolic-co-caprolactone sutures)







MONOFIT-S

Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	3-0 (2)	2-0 (3)	1-0 (3.5)	1 (4)
<b>29 m/m</b>								
1/2								
		Violet	70cm		MVS329	MVS229	MVS111	
<b>31 m/m</b>								
1/2								
		Violet	70cm			MVS207	MVS103	
<b>33 m/m</b>								
1/2								
		Undyed	75cm			MWS234		
		Violet	75cm			MVS234		
<b>36 m/m</b>								
1/2								
			70cm			MVS209	MVS105	MVS0101
		Violet	90cm				MVS104	

# MONOFIT-S

(Monofilament Polyglycolic-co-caprolactone sutures)





MONOFIT-S

Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	3-0 (2)	2-0 (3)	1-0 (3.5)	1 (4)
<b>37 m/m</b>								
3/8								
		Undyed	70cm			MWS239		
1/2								
		Undyed	70cm	MWS340	MWS240			MWS0122
			75cm				MWS122	
								
		Violet	70cm	MVS340	MVS240			MVS0122
			75cm				MVS122	
<b>41 m/m</b>								
1/2								
		Undyed	70cm			MWS244	MWS127	MWS0127
		Violet	70cm			MVS244	MVS127	MVS0127

# MONOFIT-S

(Monofilament Polyglycolic-co-caprolactone sutures)

MONOFIT-S


Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	3-0 (2)	2-0 (3)	1-0 (3.5)	1 (4)
<b>63 m/m</b>								
3/8		Violet	70cm					MVS0163
								
<b>85 m/m</b>								
3/8		Violet	70cm					MVS0185
								

## Non-neededled Suture

Suture		Product Code							
Color	Length	USP (EP)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)	1-0 (3.5)	1 (4)
Undyed Violet	30cm x 1	MWS 6301	MWS 5301	MWS 4301	MWS 3301	MWS 2301	MWS 1301	MWS 01301	
		MVS 6301	MVS 5301	MVS 4301	MVS 3301	MVS 2301	MVS 1301	MVS 01301	
Undyed Violet	45cm x 1	MWS 6451	MWS 5451	MWS 4451	MWS 3451	MWS 2451	MWS 1451	MWS 01451	
		MVS 6451	MVS 5451	MVS 4451	MVS 3451	MVS 2451	MVS 1451	MVS 01451	
Undyed Violet	75cm x 1	MWS 6751	MWS 5751	MWS 4751	MWS 3751	MWS 2751	MWS 1751	MWS 01751	
		MVS 6751	MVS 5751	MVS 4751	MVS 3751	MVS 2751	MVS 1751	MVS 01751	
Undyed Violet	100cm x 1	MWS 61001	MWS 51001	MWS 41001	MWS 31001	MWS 21001	MWS 11001	MWS 011001	
		MVS 61001	MVS 51001	MVS 41001	MVS 31001	MVS 21001	MVS 11001	MVS 011001	
Undyed Violet	150cm x 1	MWS 61501	MWS 51501	MWS 41501	MWS 31501	MWS 21501	MWS 11501	MWS 011501	
		MVS 61501	MVS 51501	MVS 41501	MVS 31501	MVS 21501	MVS 11501	MVS 011501	

# MONOFIT-L

(Monofilament  
Polydioxanone sutures)



## Monofit-L (모노피트-L 흡수성봉합사)

PDO (Polydioxanone Suture)

### DESCRIPTION 설명

촉매반응으로 나타난 Para-Dioxanone의 단위체를 중합한 형태의 흡수성봉합사로 가수분해 되어 지며, 180일에서 240일 안에 완전히 흡수된다.

한 가닥으로 구성된 단선사로서 우수한 유연성을 가지고 있고, 표면이 부드러워 조직 통과 시 낮은 자극 및 상처발생을 최소화한다. 또한 높은 초기 인장강도로 상처가 치유되는 기간 동안 강한 힘으로 지탱해주고 그 후에는 빠르게 흡수된다.

Polymer absorbable suture, which is made from Para-Dioxanone monomer through catalyst reaction, can be hydrolyzed and fully absorbable within 180 ~240 days and it has superior flexibility with its single strand suture and its smooth surface ensures less trauma in the tissue by easy passage-through and minimum injury occurrence. And, its very high initial tensile strength preserves the stitched part strongly during healing period, and later it shows fast absorption.

### COLOUR 색상

Violet, Beige(Undyed)

### STERILIZATION 멸균

Ethylene Oxide

### SIZE 사이즈

USP 7/0 ~ 2

EP 0.5 ~ 5

% of original strength remaining

2 Weeks	85%
3 Weeks	80%
4 Weeks	70%








### CERTIFICATION

- Certificate of GMP (KFDA)
- ISO9001
- ISO13485

※Catalogue에 명시된 제품 외에도 주문에 따라 생산 가능 합니다.

# MONOFIT-L













(Monofilament Polydioxanone sutures)

Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	7-0 (0.5)	6-0 (0.7)	5-0 (1)	4-0 (1.5)
<b>8m/m</b>								
3/8 ☉		Violet	25cm	MVL702				
			45cm		MVL601			
3/8 ☉ Loop		Violet	150cm		MVL602L			
3/8 ☒		Violet	25cm	MVL701				
			45cm				MVL500	
<b>9m/m</b>								
3/8 ☉		Violet	75cm			MVL603		
3/8 ☉		Violet	70cm	MVL703TDN				
3/8 ☒		Violet	45cm			MVL6110	MVL510	
1/2 ☉		Violet	30cm				MVL511	

MONOFIT-L

# MONOFIT-L













(Monofilament Polydioxanone sutures)

Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	7-0 (0.5)	6-0 (0.7)	5-0 (1)	4-0 (1.5)
<b>11 m/m</b>								
								
		Violet	45cm					MVL4131DN
3/8								
		Violet	45cm			MVL605		
			70cm			MVL606		
3/8								
		Undyed	45cm			MWL607		
		Violet	45cm			MVL607		
3/8								
		Undyed	45cm			MWL6071	MWL517	
		Violet	45cm			MVL617	MVL517 MVL517P	
3/8								
		Violet	45cm			MVL617DN		
1/2								
		Violet	45cm			MVL613		
			30cm				MVL513	



# MONOFIT-L

(Monofilament Polydioxanone sutures)













Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	7-0 (0.5)	6-0 (0.7)	5-0 (1)	4-0 (1.5)
<b>13m/m</b>								
3/8 		Undyed	70cm				MWL502	
		Violet	45cm			MVL610		
			60cm					MVL421R
			70cm	MVL704			MVL502	
3/8 		Undyed	45cm			MWL6081	MWL5031	
		Violet	45cm			MVL621	MVL5031	MVL421
3/8 		Undyed	45cm			MWL6082	MWL5032	
		Violet	45cm			MVL6082	MVL5032	
3/8 		Undyed	45cm			MWL608	MWL503	
		Violet	45cm			MVL608	MVL503	
1/2 		Violet	45cm			MVL609	MVL501	
1/2 		Undyed	45cm				MWL519	
		Violet	45cm			MVL619		

MONOFIT-L

# MONOFIT-L









(Monofilament Polydioxanone sutures)

MONOFIT-L

Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)
<b>14m/m</b>								
1/2		Violet	45cm				MVL401	
								
<b>15m/m</b>								
3/8		Violet	70cm				MVL426R	
							MVL426PR	
3/8		Violet	45cm			MVL526	MVL426	
								
3/8		Violet	70cm				MVL426DN	
								
<b>17m/m</b>								
3/8		Undyed	70cm				MWL404	
		Violet	70cm				MVL507	
3/8		Undyed	45cm				MWL506	MWL403
			70cm				MWL505	MWL402
			45cm				MVL506	MVL403
			70cm				MVL505	MVL402
			Violet	45cm	MVL506	MVL403		
			70cm	MVL505	MVL402			

# MONOFIT-L

(Monofilament Polydioxanone sutures)











Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)
<b>17m/m</b>								
3/8 		Undyed	45cm				MWL4031	
		Violet	45cm			MVL5061		
1/2 		Violet	45cm			MVL529	MVL429	
1/2 		Violet	90cm			MVL508DN	MVL405DN	
<b>18m/m</b>								
3/8 		Undyed	45cm				MWL434	
		Violet	45cm				MVL434	

MONOFIT-L

# MONOFIT-L





(Monofilament Polydioxanone sutures)

MONOFIT-L

Needle		Suture		Product Code					
Length	Description	Color	Length	USP (EP)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)
<b>18m/m</b>									
3/8			45cm				MVL406		
		Violet	70cm				MVL407		
1/2			45cm	MVL630	MVL530	MVL430	MVL308		
		Violet	45cm						
<b>19m/m</b>									
1/2			45cm		MWL535	MWL435			
		Undyed	45cm						
<b>21 m/m</b>									
3/8			45cm				MVL441	MVL319	
		Violet	45cm						
1/2			70cm				MWL442	MWL320	MWL220
		Undyed	70cm						
		Violet	70cm				MVL408	MVL301	

# MONOFIT-L

(Monofilament Polydioxanone sutures)







Needle		Suture		Product Code							
Length	Description	Color	Length	USP (EP)	4-0 (1.5)	3-0 (2)	2-0 (3)	1-0 (3.5)	1 (4)		
<b>25 m/m</b>											
3/8 		Undyed	45cm			MWL303					
			70cm			MWL324					
		Violet	45cm	MVL409	MVL303						
			70cm		MVL324	MVL224					
		1/2 		Violet	70cm	MVL447	MVL304	MVL202	MVL1071		
				Violet	45cm		MVL302				
70cm			MVL201								

MONOFIT-L

# MONOFIT-L





(Monofilament Polydioxanone sutures)

MONOFIT-L

Needle		Suture		Product Code					
Length	Description	Color	Length	USP (EP)	4-0 (1.5)	3-0 (2)	2-0 (3)	1-0 (3.5)	1 (4)
<b>28 m/m</b>									
1/2									
		Violet	70cm		MVL410				MVL108
<b>31 m/m</b>									
1/2									
		Undyed	70cm						MWL101
		Violet	70cm				MVL204	MVL101	
<b>37 m/m</b>									
1/2									
		Violet	70cm				MVL240	MVL122	MVL0122

# MONOFIT-L

(Monofilament Polydioxanone sutures)





Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	2-0 (3)	1-0 (3.5)	1 (4)	2 (5)
<b>41 m/m</b>								
1/2		Undyed	70cm			MWL106		
			70cm			MVL106	MVL0127	MVL0227
		Violet	90cm				MVL0103	MVL0201
1/2	 Loop	Violet	150cm			MVL105L	MVL0104L	
1/2		Violet	90cm				MVL0102	MVL0202

MONOFIT-L

# MONOFIT-L

(Monofilament Polydioxanone sutures)


MONOFIT-L

Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	2-0 (3)	1-0 (3.5)	1 (4)	2 (5)
<b>49 m/m</b>								
1/2		Violet	90cm				MVL0106	MVL0203
								
1/2	 Loop	Violet	150cm			MVL107L	MVL0107L	
								

## Non-needed Suture

Suture		Product Code									
Color	Length	USP (EP)	7-0 (0,5)	6-0 (0,7)	5-0 (1)	4-0 (1,5)	3-0 (2)	2-0 (3)	1-0 (3,5)	1 (4)	2 (5)
Violet	30cm x 1	MVL 7301	MVL 6301	MVL 5301	MVL 4301	MVL 3301	MVL 2301	MVL 1301	MVL 01301	MVL 02301	
Violet	45cm x 1	MVL 7451	MVL 6451	MVL 5451	MVL 4451	MVL 3451	MVL 2451	MVL 1451	MVL 01451	MVL 02451	
Violet	75cm x 1	MVL 7751	MVL 6751	MVL 5751	MVL 4751	MVL 3751	MVL 2751	MVL 1751	MVL 01751	MVL 02751	
Violet	100cm x 1	MVL 71001	MVL 61001	MVL 51001	MVL 41001	MVL 31001	MVL 21001	MVL 11001	MVL 011001	MVL 021001	
Violet	150cm x 1	MVL 71501	MVL 61501	MVL 51501	MVL 41501	MVL 31501	MVL 21501	MVL 11501	MVL 011501	MVL 021501	





# CHROMIC

(Catgut)

CHROMIC

## Chromic Catgut (크로믹 장선 흡수성봉합사)

Collagen Suture

### DESCRIPTION 설명

소의 Serosa Layer나 양의 Submucosal Layer에서 추출한 콜라겐으로 만들어진 흡수성 봉합사로 효소반응에 의해 분해되며, 흡수되는데는 약 90일 정도 소요된다.  
한 가닥으로 구성된 단선사로서 표면이 부드러워 조직을 부드럽게 통과하며 매듭의 안전성이 뛰어나다. 또한 강도를 높이기 위하여 Tru-Chromicizing 처리를 하였고, 약 3주에서 4주 정도 상처를 지지해 줄 수 있다.

Collagen absorbable suture, which is made of purified beef or sheep intestines (Serosa Layer or submucosal layer), can be breakdown through enzyme reaction and fully absorbable within 90 days and it has superior flexibility with its single strand suture, it's smooth surface ensures less trauma in the tissue by easy passage-through and superior secure knotting. And, it is treated with Tru-chromicizing to enhance tensile strength, so it ensures 3~4 weeks strength retention.

### COLOUR 색상

Brown

### STERILIZATION 멸균

Irradiation

### SIZE 사이즈

USP 6/0 ~ 3

EP 1 ~ 7





















### CERTIFICATION

- Certificate of GMP (KFDA)
- ISO9001
- ISO13485

※Catalogue에 명시된 제품 외에도 주문에 따라 생산 가능 합니다.

# CHROMIC

(Catgut)

Needle		Suture	Product Code					
Length	Description	Length	USP (EP)	6-0 (1)	5-0 (1.5)	4-0 (2)	3-0 (3)	2-0 (3.5)
<b>8m/m</b>								
1/4 		50cm				C480DN		
3/8 		50cm	C607					
<b>9m/m</b>								
3/8 		50cm	C610	C510				
1/2 		30cm	C611	C511				
<b>11m/m</b>								
		50cm				C4131DN		
3/8 		50cm	C617	C517	C417			
1/2 		30cm	C613	C513	C413			
<b>13m/m</b>								
3/8 		60cm			C421R			
3/8 		50cm	C621	C521	C421			
1/2 		50cm		C519				











# CHROMIC

(Catgut)

Needle		Suture	Product Code					
Length	Description	Length	USP (EP)	6-0 (1)	5-0 (1.5)	4-0 (2)	3-0 (3)	2-0 (3.5)
<b>14m/m</b>								
1/2	⊙	30cm	C622	C522				
	⌒	50cm				C422		
<b>15m/m</b>								
3/8	⊙	75cm				C426R		
3/8	▽	50cm	C626	C526	C426	C304		
<b>17m/m</b>								
3/8	▽	50cm			C528			
	⌒	75cm				C428	C306	
1/2	⊙	50cm			C529	C429	C307	
<b>18m/m</b>								
3/8	▽	50cm				C434	C312	
1/2	⊙	50cm	C630	C530	C430	C308		











# CHROMIC

(Catgut)

Needle		Suture	Product Code					
Length	Description	Length	USP (EP)	5-0 (1.5)	4-0 (2)	3-0 (3)	2-0 (3.5)	1-0 (4)
<b>19m/m</b>								
1/2		50cm	C535	C435	C313			
								
<b>21m/m</b>								
3/8		50cm		C441	C319			
								
1/2		50cm	C542					
		75cm		C442	C320 C320P	C220	C120	
1/2		50cm		C440	C318			
								
<b>23m/m</b>								
1/2		75cm		C443	C321	C221		
								











# CHROMIC

(Catgut)

Needle		Suture	Product Code					
Length	Description	Length	USP (EP)	4-0 (2)	3-0 (3)	2-0 (3.5)	1-0 (4)	1 (5)
<b>25 m/m</b>								
3/8 		75cm	C4465		C324	C224		
1/2 		75cm	C447		C325 C325P	C225 C225P	C107	
<b>28 m/m</b>								
1/2 		75cm			C326	C226		
<b>29 m/m</b>								
1/2 		75cm			C329	C229	C111	
<b>31 m/m</b>								
1/2 		75cm			C330	C230	C112	C0112









CHROMIC

# CHROMIC | (Catgut)

Needle		Suture	Product Code					
Length	Description	Lenght	USP (EP)	3-0 (3)	2-0 (3.5)	1-0 (4)	1 (5)	2 (6)
<b>33 m/m</b>								
1/2								
		75cm			C234	C116		
<b>36 m/m</b>								
1/2								
		75cm			C236	C118	C0118	
<b>37 m/m</b>								
3/8								
		75cm			C239			
1/2								
		75cm	C340	C240	C122	C0122		
<b>38 m/m</b>								
1/2								
		75cm	C341	C241				

# CHROMIC | (Catgut)

CHROMIC

Needle		Suture	Product Code					
Length	Description	Length	USP (EP)	3-0 (3)	2-0 (3.5)	1-0 (4)	1 (5)	2 (6)
<b>41 m/m</b>								
1/2		75cm			C244	C127	C0127	
						C126	C0126	
<b>45 m/m</b>								
1/2		75cm				C132		C0232
								
<b>63 m/m</b>								
3/8		75cm						C0163
								
<b>85 m/m</b>								
3/8		75cm						C0185
								

## Non-needed Suture

Suture	Product Code							
Length	USP (EP)	5-0 (1.5)	4-0 (2)	3-0 (3)	2-0 (3.5)	1-0 (4)	1 (5)	2 (6)
75cm x 1		C 57501	C 47501	C 37501	C 27501	C 17501	C 017501	C 027501
150cm x 1		C 51501	C 41501	C 31501	C 21501	C 11501	C 011501	C 021501





# SILK

(Silicone Coated Braided)

SILK

## Silk (실크 비흡수성봉합사)

### DESCRIPTION 설명

Degumming 공법을 통하여 99% 이상의 단백질로 구성되어 있는 순도가 높은 천연 실크로 만들어진 비흡수성봉합사로 약 3개월 이상 오랜 기간 상처를 지탱할 수 있다. 여러 가닥으로 구성된 복선사로서 표면이 부드러워 조직을 부드럽게 통과하며 매듭의 안전성이 뛰어나다. 또한 실리콘 코팅 처리되어 조직을 부드럽게 통과하고 복선사의 단점인 모세관현상을 줄여 박테리아에 감염될 위험을 최소화 한다. 일반적으로 피부봉합이나 결찰에 많이 이용되고 그 외에도 다양한 용도로 사용된다.

Non-absorbable suture, which is made of high purity natural silk, consisted with 99% protein through degumming method, it ensures about 3 months strength retention in tissue. It is made from multi strand and coated with silicone, it ensures the smooth surface and less trauma in the tissue by easy passage-through and superior knot stability. Also it lessens the infection risk by minimizing the capillary phenomenon that is the weak point of multi stranded suture. It is widely used on general skin suturing works, deligation and various propose.

### COLOUR 색상

Black, White(Undyed)

### STERILIZATION 멸균

Ethylene Oxide

### SIZE 사이즈










USP 8/0 ~ 2











EP 0.4 ~ 5









### CERTIFICATION


- Certificate of GMP (KFDA)
- ISO13485
- ISO9001
- CE Mark







※Catalogue에 명시된 제품 외에도 주문에 따라 생산 가능 합니다.



Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	8-0 (0.4)	7-0 (0.5)	6-0 (0.7)	5-0 (1)
<b>5m/m</b>								
3/8 		Black	30cm		SK818			
<b>6m/m</b>								
3/8 		Black	50cm			SK721		
3/8 		White	50cm			SW722		
		Black	50cm			SK722		
		Black	30cm			SK7221		
	 220 Microns				SK822			
<b>7m/m</b>								
3/8 		Black	30cm		SK823			
			50cm			SK723		
3/8 		Black	30cm		SK824			
			50cm			SK724		
3/8 		Black	50cm			SK723DN		
3/8 		Black	50cm			SK724DN		

Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	8-0 (0.4)	7-0 (0.5)	6-0 (0.7)	5-0 (1)
<b>8m/m</b>								
1/4 		Black	50cm				SK6081DN	
3/8 		Black	50cm	SK8251	SK725		SK6071 SK607	
3/8 	 220 Microns	Black	50cm			SK7261		
3/8 		Black	50cm	SK826		SK726 SK726P	SK608	
3/8 	 220 Microns	Black	50cm			SK7261DN		
3/8 		Black	50cm			SK726DN	SK608DN	






Needle		Suture		Product Code				
Length	Description	Color	Lenght	USP (EP)	8-0 (0.4)	7-0 (0.5)	6-0 (0.7)	5-0 (1)
<b>9m/m</b>								
3/8 			50cm			SK727	SK609	
		Black	75cm					SK5171
3/8 		White	50cm			SW728	SW610	
		Black	50cm			SK728	SK610	
3/8 		Black	50cm			SK728DN	SK610DN	
1/2 		Black	30cm			SK729	SK611	









Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	7-0 (0.5)	6-0 (0.7)	5-0 (1)	4-0 (1.5)
<b>10.5m/m</b>								
3/8		Black	75cm				SKJ501	
								
<b>11 m/m</b>								
3/8		White	50cm			SW617	SW517	
		Black	50cm	SK717		SK617 SK617P	SK517 SK517P	SK417
			75cm			SK6175		
3/8		Black	50cm			SK617RDN	SK517RDN	
3/8		Black	50cm			SK617DN SK617PDN	SK517DN	
1/2		Black	30cm	SK731		SK613	SK513	SK413
1/2		Black	50cm			SK613DN		
1/2		Black	50cm	SK7161DN				









Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	7-0 (0.5)	6-0 (0.7)	5-0 (1)	4-0 (1.5)
<b>13m/m</b>								
3/8		Black	50cm			SK621R	SK521R	SK421R
			75cm				SK5215R	SK421PR
3/8		Black	50cm			SK621	SK521	SK421
			75cm				SK5211	
3/8		Black	50cm					SK421DN
<b>14m/m</b>								
1/2		Black	30cm			SK622	SK522	SK422











Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)
<b>15m/m</b>								
3/8	⊙	Black	50cm			SK526R		
			75cm				SK426R SK426PR	
3/8	▽	Black	50cm				SK425	
3/8	▽	Black	50cm	<div style="border: 1px dashed black; padding: 2px;">                     White 45cmx1 Black 45cmx3                 </div>			SKW426M	
					SK626 SK626P	SK526 SK526P	SK426 SK426P	SK304
<b>16m/m</b>								
3/8	⊙	Black	75cm				SKJ401	



















Needle		Suture		Product Code				
Length	Description	Color	Lenght	USP (EP)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)
<b>17m/m</b>								
3/8		Black	50cm	SK528	SK428	SK306		
			75cm				SK206	
1/2		Black	45cm x 7		SK4181			
			50cm	SK529	SK429	SK307		
			75cm				SK307P	
1/2		Black	50cm	SK527				







Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)
<b>18m/m</b>								
3/8 		White	50cm			SW434	SW312	
		Black	50cm	SK534		SK434 SK434P	SK312	
1/2 			45cm x 5			SK430RM		
		Black	50cm			SK430	SK308	SK213
<b>19m/m</b>								
3/8 			50cm	SK5351			SK3131	
		Black	75cm			SK4351		
1/2 		Black	50cm	SK535		SK435	SK313	

Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)
<b>21 m/m</b>								
3/8								
		Black	75cm				SK3201	
3/8								
		Black	50cm	SK539P		SK441 SK439	SK319 SK317	SK219 SK217
1/2								
		Black	45cm x 5				SK320RM	
			75cm	SK542		SK442	SK320	SK220
1/2								
		Black	50cm				SK318	SK218

Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	4-0 (1.5)	3-0 (2)	2-0 (3)	1-0 (3.5)
<b>23m/m</b>								
1/2								
		Black	75cm	SK443		SK321		
<b>24m/m</b>								
1/2								
		Black	50cm			SK322		
<b>25m/m</b>								
3/8								
		White	75cm			SW324		
		Black	50cm	SK446		SK324 SK324P	SK224 SK224P	
1/2								
			45cm x 8	SK447RM		SK325RM	SK225RM	
		Black	75cm x 8	SK4471RM		SK3251RM	SK2251RM	
			75cm	SK447		SK325 SK325P SKJ301	SK225 SK225P	SK107
1/2								
		Black	75cm				SK223	

Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	3-0 (2)	2-0 (3)	1-0 (3.5)	1 (4)
<b>28 m/m</b>								
1/2								
		Black	75cm		SK326			
<b>29 m/m</b>								
3/8								
		White	75cm			SW228	SW110	
		Black	75cm		SK328	SK228	SK110	SK0110
1/2								
		Black	75cm		SK329	SK229	SK111	
<b>31 m/m</b>								
1/2								
		Black	75cm				SK112	SK0112


Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	3-0 (2)	2-0 (3)	1-0 (3.5)	1 (4)
<b>33 m/m</b>								
3/8								
		White	75cm					SW0115
		Black	75cm			SK233	SK115	SK0115
1/2								
		Black	75cm	SK334		SK234	SK116	SK0116
<b>37 m/m</b>								
3/8								
		Black	75cm			SK239	SK121	SK0121
1/2								
		Black	75cm			SK240	SK122	SK0122

Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	2-0 (3)	1-0 (3.5)	1 (4)	2 (5)
<b>41 m/m</b>								
3/8								
		Black	75cm			SK125	SK0125	
			100cm		SK243			
1/2								
		Black	75cm		SK244 SK244P	SK126	SK0126	SK0227
<b>45 m/m</b>								
3/8								
		Black	75cm				SK0131	

## Non-needed Suture

Suture		Product Code									
Color	Length	USP (EP)	7-0 (0.5)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)	1-0 (3.5)	1 (4)	2 (5)
White Black	40cm x 4					SW44004 SK 44004	SW34004 SK 34004				
White Black	45cm x 7									SW01457 SK 01457	
White Black	45cm x 10			SW64510 SK 64510	SW54510 SK 54510	SW44510 SK 44510	SW34510 SK 34510	SW24510 SK 24510	SW14510 SK 14510	SW014510 SK 014510	SW024510 SK 024510
White Black	45cm x 17			SW64517 SK 64517	SW54517 SK 54517	SW44517 SK 44517	SW34517 SK 34517	SW24517 SK 24517	SW14517 SK 14517	SW014517 SK 014517	SW024517 SK 024517
White Black	50cm x 10	WK75010 SK 75010									
White Black	60cm x 13			SW66013 SK 66013	SW56013 SK 56013	SW46013 SK 46013	SW36013 SK 36013	SW26013 SK 26013	SW16013 SK 16013	SW016013 SK 016013	SW026013 SK 026013
White Black	75cm x 10	SW77510 SK 77510		SW67510 SK 67510	SW57510 SK 57510	SW47510 SK 47510	SW37510 SK 37510	SW27510 SK 27510	SW17510 SK 17510	SW017510 SK 017510	SW027510 SK 027510
White Black	110cm x 7			SW61107 SK 61107	SW51107 SK 51107	SW41107 SK 41107	SW31107 SK 31107	SW21107 SK 21107	SW11107 SK 11107	SW011107 SK 011107	SW021107 SK 021107
White Black	150cm x 1			SW61501 SK 61501	SW51501 SK 51501	SW41501 SK 41501	SW31501 SK 31501	SW21501 SK 21501	SW11501 SK 11501	SW011501 SK 011501	SW021501 SK 021501
White Black	180cm x 1					SW41801 SK 41801	SW31801 SK 31801	SW21801 SK 21801	SW11801 SK 11801	SW011801 SK 011801	SW021801 SK 021801





# NYLON

(Monofilament)

NYLON

## Nylon (나일론 비흡수성봉합사)

Polyamide

### DESCRIPTION 설명

Polyamide 6과 Polyamide 6.6으로 만들어진 비흡수성봉합사로, 조직내에서 영구적인 인장강도를 유지한다. 단선사로 표면이 매끄러워 조직 투과력이 뛰어나고 박테리아의 증식을 최소화 할 수 있어 감염의 위험이 적다. 또한 단선사이지만 매우 유연하고 핸들링이 뛰어나 주로 피부봉합에 많이 사용되며, 재질의 특성상 매우 가는 봉합사로도 생산 가능하여 안과 수술이나 미세 수술에도 많이 사용된다.

Non-absorbable suture, which is made of Polyamide 6 and Polyamide 6.6, ensures prolonged permanent tensile strength retention in tissue. It is made from single strand, so it ensures the smooth surface and less trauma in the tissue by easy passage-through, and it lessens the infection risk by minimizing bacterial growth. Also, it's flexibility provides ensures good handling for skin suture, the characteristics of nylon material makes it possible to manufacture very fine thread, so it is widely used for ophthalmic surgery and microsurgery.

### COLOUR 색상

Black, Blue, White(Undyed)

### STERILIZATION 멸균

Ethylene Oxide

### SIZE 사이즈

USP 11/0 ~ 2

EP 0.1 ~ 5

### CERTIFICATION

- Certificate of GMP (KFDA)
- ISO13485
- ISO9001
- CE Mark

※Catalogue에 명시된 제품 외에도 주문에 따라 생산 가능 합니다.

# NYLON

(Monofilament)

Needle		Suture		Product Code					
Length	Description	Color	Lenght	USP (EP)	11-0 (0.1)	10-0 (0.2)	9-0 (0.3)	8-0 (0.4)	7-0 (0.5)
<b>3m/m</b>									
3/8 ☉	75 Microns	Black	15cm	NK1102	NK1002				
<b>3.5m/m</b>									
3/8 ☉	75 Microns	Black	15cm	NK1104	NK1004				
<b>3.75m/m</b>									
3/8 ☉	75 Microns	Black	15cm		NK1041				
<b>4m/m</b>									
3/8 ☉	100 Microns	Black	15cm		NK1007	NK907			
3/8 ☹	100 Microns	Black	15cm		NK1008	NK908			
3/8 ☹	130 Microns	Black	15cm					NK808	
3/8 ☉	150 Microns	Black	15cm					NK807	
<b>4.75m/m</b>									
3/8 ☉	130 Microns	Black	15cm		NK1031	NK9151	NK8151		








# NYLON

(Monofilament)

Needle		Suture		Product Code					
Length	Description	Color	Length	USP (EP)	11-0 (0.1)	10-0 (0.2)	9-0 (0.3)	8-0 (0.4)	7-0 (0.5)
<b>5m/m</b>									
3/8	  100 Microns	Black	15cm			NK1051	NK951		
3/8	  150 Microns	Black	15cm				NK915	NK815	
3/8	 	Blue	30cm					NB817	
3/8	  100 Microns	Black	45cm				NK951DN		
<b>6m/m</b>									
3/8	  220 Microns	Black	30cm			NK1070	NK970		
3/8	 	Blue	30cm					NB821	NB721
		Black	30cm					NK821	NK721
3/8	 	White	50cm						NW722
		Blue	30cm					NB822	
			50cm						NB722
		Black	30cm					NK822	
			50cm						NK722

# NYLON

(Monofilament)

Needle		Suture		Product Code					
Length	Description	Color	Length	USP (EP)	10-0 (0.2)	9-0 (0.3)	8-0 (0.4)	7-0 (0.5)	6-0 (0.7)
<b>6m/m</b>									
3/8		Black	30cm				NK821DN		
3/8	 150 Microns	Black	30cm	NK1090DN	NK990DN				
3/8	 220 Microns	Black	30cm	NK1070DN	NK970DN				
<b>6.5m/m</b>									
3/8	 130 Microns	Black	15cm				NK890		
<b>7m/m</b>									
3/8		Blue	50cm				NB823	NB723	
3/8		White	50cm				NW824	NW724	
		Blue	50cm				NB824	NB724	
		Black	50cm					NK724	
<b>8m/m</b>									
3/8		Blue	50cm				NB825	NB725	NB607
		Black	50cm				NK825	NK725	NK607









# NYLON

(Monofilament)

Needle		Suture		Product Code					
Length	Description	Color	Length	USP (EP)	8-0 (0.4)	7-0 (0.5)	6-0 (0.7)	5-0 (1)	4-0 (1.5)
<b>8 m/m</b>									
3/8		White	50cm			NW726 NW726P	NW608 NW608P	NW508 NW508P	
		Blue	50cm	NB826		NB726 NB726P	NB608 NB608P		
		Black	50cm	NK826		NK726 NK726P	NK608		
3/8		Black	50cm			NK7261			
3/8		Black	50cm	NK825PDN			NK607DN		
3/8		Blue	50cm			NB726DN			
1/2	220 Microns	Black	50cm	NK826P					
<b>9 m/m</b>									
3/8		White	50cm			NW727			
		Blue	50cm			NB727	NB609		
3/8		White	50cm			NW728	NW610		
		Blue	50cm			NB728	NB610		
		Black	50cm			NK728	NK610		
3/8		Blue	50cm				NB6101P		
		Black	50cm				NK6101P	NK5101P	













# NYLON

(Monofilament)

Needle		Suture		Product Code					
Length	Description	Color	Length	USP (EP)	8-0 (0.4)	7-0 (0.5)	6-0 (0.7)	5-0 (1)	4-0 (1.5)
<b>9 m/m</b>									
3/8 ☒		White	50cm			NW728DN			
		Blue	50cm			NB728DN	NB610DN		
		Black	50cm			NK728DN	NK610DN		
1/2 ⊙		Blue	30cm	NB829	NB729	NB611			
<b>11 m/m</b>									
3/8 ⊙		White	50cm			NW717R	NW617R	NW517R	
		Blue	50cm			NB717R	NB617R	NB517R	
3/8 ☒		White	50cm			NW717 NW717P	NW617 NW617P	NW517 NW517P	NW417
		Blue	50cm			NB717 NB717P	NB617 NB617P	NB517 NB517P	NB417 NB417P
		Black	50cm			NK717 NK717P	NK617 NK617P	NK517 NK517P	NK417P
3/8 ☒		Blue	50cm				NB616		
3/8 ☒		White	50cm					NW516DN	
3/8 ☒		Blue	50cm			NB717DN	NB617DN		
		Black	50cm			NK717DN			
1/2 ⊙		Blue	30cm				NB613	NB513	NB413

# NYLON










(Monofilament)

Needle		Suture		Product Code						
Length	Description	Color	Length	USP (EP)	8-0 (0.4)	7-0 (0.5)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	
<b>11 m/m</b>										
1/2		White	50cm					NW514		
										
1/2		White	50cm				NW615 NW615P	NW515P		
		Blue	50cm				NB615			
		Black	50cm				NK615			
<b>13 m/m</b>										
3/8		White	50cm					NW621 NW621P	NW521 NW521P	NW421 NW421P
										
		Black	50cm				NK621	NK521	NK421	
3/8		White	50cm			NW621DN				
		Blue	50cm			NB621DN				
		Black	50cm			NK621DN				
3/8		Black	45cm		NK720D					
1/2		White	50cm				NW619P	NW519P		
										
		Blue	50cm				NB619	NB519		
1/2		Blue	50cm				NB518			















# NYLON

(Monofilament)

Needle		Suture		Product Code					
Length	Description	Color	Length	USP (EP)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)
<b>14m/m</b>									
1/2		Blue	30cm	NB622	NB522	NB422	NB300		
									
<b>15m/m</b>									
3/8		Blue	50cm		NB526R				
3/8		White	50cm	NW626	NW526 NW526P	NW426	NW304P		
		Blue	50cm	NB626 NB626P	NB526 NB526P	NB426	NB304		
		Black	35cm	NK626					
50cm	NK626P		NK526 NK526P	NK426					
1/2		Blue	50cm		NB524	NB424			
<b>17m/m</b>									
3/8		Blue	50cm		NB528 NB528P	NB428	NB306	NB206	
		Black	50cm		NK528				
1/2		Blue	35cm			NB429L			
1/2		Blue	50cm		NB529	NB429	NB307	NB207	
1/2		Blue	50cm		NB527	NB427			







# NYLON

(Monofilament)

Needle		Suture		Product Code							
Length	Description	Color	Length	USP (EP)	8-0 (0.4)	7-0 (0.5)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)
<b>18m/m</b>											
3/8		White	50cm						NW434	NW312	
		Blue	50cm				NB634	NB534	NB434 NB434P	NB312 NB312P	
		Black	50cm						NK434	NK312	
3/8		White	50cm				NW634DN				
1/2		Blue	50cm				NB630	NB530	NB430	NB308	
<b>19m/m</b>											
3/8		White	45cm			NW7019DN					
		Blue	45cm	NB8019WR	NB7019DN						
		Black	45cm		NK7019DN						
1/2		White	50cm			NW735P					
		Blue	50cm					NB535	NB435	NB313	
<b>21 m/m</b>											
3/8		White	50cm						NW439	NW317	
		Blue	50cm					NB541	NB441 NB439	NB319 NB317	NB219 NB217
		Black	50cm							NK317	
1/2		Blue	50cm					NB542	NB442	NB320	NB220
1/2		Blue	50cm					NB540	NB440	NB318	NB218















# NYLON

(Monofilament)

Needle		Suture		Product Code							
Length	Description	Color	Length	USP (EP)	8-0 (0.4)	7-0 (0.5)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)
<b>23 m/m</b>											
1/2		Blue	50cm					NB543			
			75cm						NB443	NB321	
<b>24 m/m</b>											
3/8		White	45cm		NW7024DN						
		Blue	45cm		NB7024DN						
		Black	45cm		NK7024DN						
1/2		Blue	50cm							NB322	
<b>25 m/m</b>											
3/8		White	75cm							NW324	NW224
		Blue	50cm					NB546	NB446	NB3245	NB324P
			75cm							NB324	NB224
		Black	75cm							NK324	NK224
3/8		Blue	75cm							NB324DN	
1/2		Blue	75cm						NB447	NB325	NB225









# NYLON

(Monofilament)

Needle		Suture		Product Code						
Length	Description	Color	Length	USP (EP)	4-0 (1.5)	3-0 (2)	2-0 (3)	1-0 (3.5)	1 (4)	2 (5)
<b>29 m/m</b>										
3/8		Blue	75cm			NB328	NB228	NB110	NB0110	NB0210
			75cm				NB229	NB111		
1/2		Blue	75cm				NB229	NB111		
			100cm				NB2291			
<b>31 m/m</b>										
1/2		Blue	75cm				NB231			
										
<b>33 m/m</b>										
3/8		Blue	75cm			NB333	NB233	NB115	NB0115	
										
1/2		Blue	75cm				NB234			
										
<b>37 m/m</b>										
3/8		Blue	75cm			NB339	NB239	NB121	NB0121	
		Black	75cm				NK239			
1/2		Blue	75cm					NB122	NB0122	
										

# NYLON











(Monofilament)

Needle		Suture		Product Code						
Length	Description	Color	Length	USP (EP)	4-0 (1.5)	3-0 (2)	2-0 (3)	1-0 (3.5)	1 (4)	2 (5)
<b>41 m/m</b>										
3/8 		Blue	75cm				NB243	NB125	NB0125	NB0225
1/2 		Blue	75cm					NB126	NB0126	
1/2 		Blue	75cm					NB124	NB0124	
		Black	75cm					NK124		
<b>45 m/m</b>										
3/8 			75cm					NB131	NB0131	NB0231
		Blue	100cm				NB249			

NYLON







# NYLON

(Monofilament)

Needle		Suture		Product Code					
Length	Description	Color	Length	USP (EP)	3-0 (2)	2-0 (3)	1-0 (3.5)	1 (4)	2 (5)
<b>46 m/m</b>									
1/2		Blue	75cm						NB0233
									
<b>49 m/m</b>									
1/2		Blue	100cm						NB0150R
									
1/2		Blue	100cm						NB0150
									
<b>55 m/m</b>									
		Blue	100cm			NB255	NB155		
									
<b>60 m/m</b>									
		Blue	75cm		NB360	NB260			
									

# NYLON


(Monofilament)

Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	2-0 (3)	1-0 (3,5)	1 (4)	2 (5)
<b>60m/m</b>								
3/8 		Blue	100cm		NB261		NB0161	
<b>93m/m</b>								
3/8 	each with 40mm silicone tubing 2 Sutuers per pack	Blue	50cm X 2				NB0193M	NB0293M
			100cm		NB193	NB0193	NB0293	
3/8 		Blue	75cm				NB0193DN	

## Non-needed Suture

Suture		Product Code									
Color	Length	USP (EP)	7-0 (0.5)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)	1-0 (3.5)	1 (4)	2 (5)
White Blue Black	40cm x 4					NW44004 NB 44004 NK 44004	NW34004 NB 34004 NK 34004				
White Blue Black	45cm x 7									NB 01457	
White Blue Black	45cm x 10			NW64510 NB 64510 NK 64510	NW54510 NB 54510 NK 54510	NW44510 NB 44510 NK 44510	NW34510 NB 34510 NK 34510	NW24510 NB 24510 NK 24510	NW14510 NB 14510 NK 14510	NB 014510	NB 024510
White Blue Black	45cm x 17			NW64517 NB 64517 NK 64517	NW54517 NB 54517 NK 54517	NW44517 NB 44517 NK 44517	NW34517 NB 34517 NK 34517	NW24517 NB 24517 NK 24517	NW14517 NB 14517 NK 14517	NB 014517	NB 024517
White Blue Black	50cm x 10	NW75010 NB 75010 NK 75010									
White Blue Black	60cm x 13			NW66013 NB 66013 NK 66013	NW56013 NB 56013 NK 56013	NW46013 NB 46013 NK 46013	NW36013 NB 36013 NK 36013	NW26013 NB 26013 NK 26013	NW16013 NB 16013 NK 16013	NB 016013	NB 026013
White Blue Black	75cm x 10			NW67510 NB 67510 NK 67510	NW57510 NB 57510 NK 57510	NW47510 NB 47510 NK 47510	NW37510 NB 37510 NK 37510	NW27510 NB 27510 NK 27510	NW17510 NB 17510 NK 17510	NB 017510	NB 027510
White Blue Black	110cm x 7			NW61107 NB 61107 NK 61107	NW51107 NB 51107 NK 51107	NW41107 NB 41107 NK 41107	NW31107 NB 31107 NK 31107	NW21107 NB 21107 NK 21107	NW11107 NB 11107 NK 11107	NB 011107	NB 021107
White Blue Black	150cm x 1			NW61501 NB 61501 NK 61501	NW51501 NB 51501 NK 51501	NW41501 NB 41501 NK 41501	NW31501 NB 31501 NK 31501	NW21501 NB 21501 NK 21501	NW11501 NB 11501 NK 11501	NB 011501	NB 021501
White Blue Black	180cm x 1					NW41801 NB 41801 NK 41801	NW31801 NB 31801 NK 31801	NW21801 NB 21801 NK 21801	NW11801 NB 11801 NK 11801	NB 011801	NB 021801





# POLYESTER

(Silicone Coated Braided)

POLYESTER

## Polyester (폴리에스터 비흡수성봉합사)

### DESCRIPTION 설명

polyester재질로 만들어진 비흡수성봉합사로, 조직내에서 영구적인 인장강도를 유지한다. 복선사로서 강도가 매우 뛰어나고 핸들링 및 매듭의 안전성이 우수하고, Silicone 코팅 처리되어 조직을 부드럽게 통과하고 매듭이 용이 하다.

Non-absorbable suture, which is made of Polyester material, ensures prolonged permanent tensile strength retention in tissue. It is made from multi strand and coated with silicone, it ensures good handling, superior secure knotting and easy passage-through and easy knot as well as strong tensile strength.

### COLOUR 색상

Green, White(Undyed)

### STERILIZATION 멸균

Ethylene Oxide

### SIZE 사이즈

USP 7/0 ~ 5

EP 0.5 ~ 7













### CERTIFICATION

- Certificate of GMP (KFDA)
- ISO13485
- ISO9001

※Catalogue에 명시된 제품 외에도 주문에 따라 생산 가능 합니다.











# POLYESTER

(Silicone coated Braided)

Needle		Suture		Product Code					
Length	Description	Color	Length	USP (EP)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)
<b>11m/m</b>									
3/8		White	75cm				PEW416RDN		
		Green	75cm				PEG416RDN		
3/8		Green	75cm	PEG617DN PEG617PDN					
									
<b>13m/m</b>									
3/8		Green	75cm			PEG5221DN	PEG4221DN		
									
<b>17m/m</b>									
1/2		Green	60cm			PEG529DN		PEG307DN	
									
1/2		Green	90cm			PEG 5291TDN	PEG 4291TDN	PEG 3071 TDN	PEG 2071TDN
									
<b>18m/m</b>									
3/8		Green	50cm				PEG 434		
									









# POLYESTER

(Silicone coated Braided)

Needle		Suture		Product Code						
Length	Description	Color	Lenght	USP (EP)	4-0 (1.5)	3-0 (2)	2-0 (3)	1-0 (3.5)	1 (4)	2 (5)
<b>21m/m</b>										
1/2 		White	55cm				PEW221			
		Green	55cm				PEG221			
		White	75cm					PEW121		
		Green	75cm					PEG121		
1/2 		Green	75cm	PEG 442DN						
<b>24m/m</b>										
1/2 		White	75cm							PEW 0224
		Green	75cm							PEG 0224
<b>25m/m</b>										
1/2 		White	75cm	PEW 447	PEW 325					
		Green	75cm	PEG 447	PEG 325					
1/2 		White	90cm				PEW 225DN			
		Green	90cm				PEG 225DN			
		White	100cm		PEW 325DN					
		Green	100cm		PEG 325DN					









# POLYESTER

(Silicone coated Braided)

Needle		Suture		Product Code						
Length	Description	Color	Length	USP (EP)	4-0 (1.5)	3-0 (2)	2-0 (3)	1-0 (3.5)	1 (4)	2 (5)
<b>25 m/m</b>										
1/2		White	90cm				PEW 2251TDN	PEW 1071TDN		
										
<b>26 m/m</b>										
1/2		Green	75cm							PEG0226
<b>28 m/m</b>										
1/2		Green	75cm					PEG108		
<b>31 m/m</b>										
1/2		White	75cm					PEW112		
		Green	75cm					PEG112		
<b>32 m/m</b>										
1/2		White	75cm						PEW0132	PEW0232
		Green	75cm						PEG0132	PEG0232





# POLYESTER

(Silicone coated Braided)

Needle		Suture		Product Code				
Length	Description	Color	Lenght	USP (EP)	1-0 (3,5)	1 (4)	2 (5)	5 (7)
<b>45 m/m</b>								
1/2		Green	75cm					PEG0245T
			75cm X4					PEG0245TM
1/2		Green	75cm					PEG0245TDN
								
<b>56 m/m</b>								
1/2		Green	75cm					PEG0556T
			75cm X4					PEG0556TM
<b>93 m/m</b>								
3/8				<b>Color</b>	<b>Lenght</b>			
		Green	50cm X2			PEG0193M		

# POLYESTER

(Silicone coated Braided)

Needle		Suture		Product Code				
Length	Description	Color	Length	USP (EP)	3-0 (2)	2-0 (3)	1-0 (3.5)	1 (4)
<b>250 m/m</b>								
Loop 								
		White	20cm	PEW 3250L				
※ 뒷면에 실제 사이즈 표시		Green	20cm	PEG 3250L				
Loop 								
		White	70cm			PEW 2250DN		
※ 뒷면에 실제 사이즈 표시		Green	70cm			PEG 2250DN		

## Non-needed Suture

Suture		Product Code									
Color	Length	USP (EP)	7-0 (0.5)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)	1-0 (3.5)	1 (4)	2 (5)
Green	40cm x 4					PEG 44004	PEG 34004				
Green	45cm x 7									PEG 01457	
Green	45cm x 10			PEG 64510	PEG 54510	PEG 44510	PEG 34510	PEG 24510	PEG 14510	PEG 014510	PEG 024510
Green	45cm x 17			PEG 64517	PEG 54517	PEG 44517	PEG 34517	PEG 24517	PEG 14517	PEG 014517	PEG 024517
Green	50cm x 10	PEG 75010									
Green	60cm x 13			PEG 66013	PEG 56013	PEG 46013	PEG 36013	PEG 26013	PEG 16013	PEG 016013	PEG 026013
Green	75cm x 10			PEG 67510	PEG 57510	PEG 47510	PEG 37510	PEG 27510	PEG 17510	PEG 017510	PEG 027510
Green	110cm x 7			PEG 61107	PEG 51107	PEG 41107	PEG 31107	PEG 21107	PEG 11107	PEG 011107	PEG 021107
Green	150cm x 1			PEG 61501	PEG 51501	PEG 41501	PEG 31501	PEG 21501	PEG 11501	PEG 011501	PEG 021501
Green	180cm x 1					PEG 41801	PEG 31801	PEG 21801	PEG 11801	PEG 011801	PEG 021801

■ 250 m/m

● Loop


20cm : PEW 3250L  
PEG 3250L

■ 250 m/m

●

70cm : PEW 2250DN  
PEG 2250DN





# POLYPROPYLENE

(Monofilament)

## Polypropylene (폴리프로필렌 비흡수성봉합사)

### DESCRIPTION 설명

Polypropylene으로 만들어진 비흡수성봉합사로, 조직을 강하고 안정적으로 지지하면서 감염된 부위에서도 강도가 약해지지 않고 변함없는 영구적 인장강도를 지닌다. 단선사로 표면이 매끄러워 조직 투과력이 뛰어나고 박테리아의 증식을 최소화 할 수 있어 감염의 위험이 적다. 또한 실의 꼬임현상이 적고 매듭이 쉬우며 안정적이다.

Non-absorbable suture, which is made of Polypropylene material, ensures prolonged permanent tensile strength retention in tissue without any showing of weakness in severe infected tissue, while it sustains the tissue strongly and stably. It is made from single strand, it has smooth surface and superior tissue penetrating ability, and it lessens the infection risk by minimizing bacterial growth. Also, it show low suture twist and easy knot and excellent knot security.

### COLOUR 색상

Blue

### STERILIZATION 멸균

Ethylene Oxide

### SIZE 사이즈

USP 11/0 ~ 3

EP 0.1 ~ 6





















### CERTIFICATION

- Certificate of GMP (KFDA)
- ISO13485
- ISO9001

※Catalogue에 명시된 제품 외에도 주문에 따라 생산 가능 합니다.



















# POLYPROPYLENE

(Monofilament)

Needle		Suture	Product Code					
Length	Description	Length	USP (EP)	11-0 (0.1)	10-0 (0.2)	9-0 (0.3)	8-0 (0.4)	7-0 (0.5)
<b>3.5m/m</b>								
3/8 	 75 Microns	15cm			PP10041			
<b>3.75m/m</b>								
3/8 	 75 Microns	15cm			PP1004			
<b>4m/m</b>								
3/8 	 100 Microns	15cm				PP907		
<b>4.5m/m</b>								
3/8 	 150 Microns	15cm				PP911	PP811	
<b>4.75m/m</b>								
3/8 	 130 Microns	15cm				PP9111	PP8111	
<b>5m/m</b>								
3/8 	 150 Microns	15cm				PP916	PP816	
<b>6m/m</b>								
3/8 	 220 Microns	25cm			PP1016	PP917		
3/8 		30cm			PP1015		PP821	
<b>6.5m/m</b>								
3/8 	 130 Microns	15cm					PP823	
3/8 	 130 Microns	45cm					PP823DN	









# POLYPROPYLENE

(Monofilament)

Needle		Suture	Product Code					
Length	Description	Length	USP (EP)	8-0 (0.4)	7-0 (0.5)	6-0 (0.7)	5-0 (1)	4-0 (1.5)
<b>8 m/m</b>								
3/8	  220 Microns	60cm			PP725			
3/8	 	60cm	PP825					
3/8	  220 Microns	60cm			PP725DN	PP625DN		
<b>9 m/m</b>								
3/8	 	60cm			PP727DN	PP609DN		
1/2	 	30cm				PP611		
<b>11 m/m</b>								
3/8	 	50cm					PP516	
3/8	 	45cm				PP617P		
		50cm					PP517	
3/8	 	60cm				PP612PDN		
1/2	 	30cm			PP731		PP513	PP413

















# POLYPROPYLENE

(Monofilament)

Needle		Suture	Product Code					
Length	Description	Length	USP (EP)	7-0 (0.5)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)
<b>11 m/m</b>								
1/2		60cm				PP513TDN		
<b>13 m/m</b>								
3/8		75cm			PP6291	PP5211		
3/8		45cm			PP618P			
		50cm				PP521		
3/8		60cm	PP729DN		PP629DN	PP529DN		
		75cm			PP6291DN			
<b>13 m/m</b>								
1/2		60cm			PP622DN			
		75cm				PP522DN		
<b>14 m/m</b>								
1/2		30cm					PP422	
<b>15 m/m</b>								
	 150 Microns	25cm	USP (EP)	10-0 (0.2)				
			PP1017DN					
3/8		45cm			PP626 PP626P	PP526 PP526P	PP426 PP426P	PP304
		75cm						PP3041P















# POLYPROPYLENE

(Monofilament)

Needle		Suture	Product Code					
Length	Description	Length	USP (EP)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)
<b>15m/m</b>								
3/8	 	90cm			PP5291BDN			
1/2	 	45cm				PP4291		
1/2	 	90cm			PP5291DN	PP4291DN	PP3071DN	
<b>16m/m</b>								
1/2	 	90cm			PP5291P			
1/2	 	75cm			PP5290DN			
		90cm			PP5291PDN			
<b>17m/m</b>								
3/8	 	50cm			PP528			
1/2	 	50cm				PP429	PP307	PP207
1/2	 	90cm			PP5292TDN	PP4292TDN	PP3072TDN	PP2072TDN











# POLYPROPYLENE

(Monofilament)

Needle		Suture	Product Code					
Length	Description	Length	USP (EP)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)
<b>18m/m</b>								
3/8 		45cm			PP534 PP534P	PP434 PP434P	PP312	
1/2 		50cm	PP630					
<b>19m/m</b>								
1/2 		50cm			PP535	PP435		PP213
<b>19m/m</b>								
1/2 		90cm				PP435DN		
<b>21 m/m</b>								
3/8 		45cm				PP439P	PP317P	
		50cm					PP316	
1/2 		50cm			PP542	PP442		
		75cm					PP320	
1/2 		45cm					PP318	

# POLYPROPYLENE









(Monofilament)

Needle		Suture	Product Code					
Length	Description	Length	USP (EP)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)
<b>23m/m</b>								
1/2								
		75cm					PP321	
<b>24m/m</b>								
1/2								
		90cm				PP422DN	PP322DN	PP2221DN
1/2								
		90cm						PP2221TDN
<b>25m/m</b>								
3/8								
		45cm		PP546P	PP446		PP324 PP324P	
		50cm						PP224
1/2								
		75cm					PP325	PP225















# POLYPROPYLENE

(Monofilament)

Needle		Suture	Product Code					
Length	Description	Length	USP (EP)	4-0 (1.5)	3-0 (2)	2-0 (3)	1-0 (3.5)	1 (4)
<b>28 m/m</b>								
1/2								
		75cm				PP226		
1/2								
		90cm	PP426DN	PP326DN	PP226DN			
<b>29 m/m</b>								
3/8								
		75cm						PP0110
<b>31 m/m</b>								
1/2								
		75cm				PP230		
		100cm					PP112	






# POLYPROPYLENE

(Monofilament)

Needle		Suture	Product Code					
Length	Description	Length	USP (EP)	3-0 (2)	2-0 (3)	1-0 (3.5)	1 (4)	2 (5)
<b>36 m/m</b>								
1/2			75cm		PP236			
1/2			100cm			PP118T		
<b>37 m/m</b>								
1/2			75cm			PP122	PP0122	
<b>41 m/m</b>								
3/8			100cm		PP243			
1/2			100cm			PP127	PP0126	
1/2			100cm			PP124		


# POLYPROPYLENE

(Monofilament)

Needle		Suture	Product Code					
Length	Description	Length	USP (EP)	3-0 (2)	2-0 (3)	1-0 (3,5)	1 (4)	2 (5)
<b>49 m/m</b>								
1/2		100cm					PP0149	
<b>55 m/m</b>								
		90cm			PP255			
<b>60 m/m</b>								
		75cm			PP2601			
<b>70 m/m</b>								
		75cm			PP270DN			
<b>75 m/m</b>								
		100cm			PP275	PP175		

## Non-neededled Suture

Suture	Product Code										
Length	USP (EP)	8-0 (0.4)	7-0 (0.5)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)	1-0 (3.5)	1 (4)	2 (5)
40cm × 10	PP 84010	PP 74010	PP 64010	PP 54010	PP 44010	PP 34010	PP 24010	PP 14010	PP 014010	PP 024010	
45cm × 10	PP 84510	PP 74510	PP 64510	PP 54510	PP 44510	PP 34510	PP 24510	PP 14510	PP 014510	PP 024510	
45cm × 17	PP 84517	PP 74517	PP 64517	PP 54517	PP 44517	PP 34517	PP 24517	PP 14517	PP 014517	PP 024517	
50cm × 10	PP 85010	PP 75010	PP 65010	PP 55010	PP 45010	PP 35010	PP 25010	PP 15010	PP 015010	PP 025010	
50cm × 15	PP 85015	PP 75015	PP 65015	PP 55015	PP 45015	PP 35015	PP 25015	PP 15015	PP 015015	PP 025015	
60cm × 10	PP 86010	PP 76010	PP 66010	PP 56010	PP 46010	PP 36010	PP 26010	PP 16010	PP 016010	PP 026010	
60cm × 13	PP 86013	PP 76013	PP 66013	PP 56013	PP 46013	PP 36013	PP 26013	PP 16013	PP 016013	PP 026013	
75cm × 10	PP 87510	PP 77510	PP 67510	PP 57510	PP 47510	PP 37510	PP 27510	PP 17510	PP 017510	PP 027510	
100cm × 1	PP 81001	PP 71001	PP 61001	PP 51001	PP 41001	PP 31001	PP 21001	PP 11001	PP 011001	PP 021001	
120cm × 1	PP 81201	PP 71201	PP 61201	PP 51201	PP 41201	PP 31201	PP 21201	PP 11201	PP 011201	PP 021201	
150cm × 1	PP 81501	PP 71501	PP 61501	PP 51501	PP 41501	PP 31501	PP 21501	PP 11501	PP 011501	PP 021501	
180cm × 1	PP 81801	PP 71801	PP 61801	PP 51801	PP 41801	PP 31801	PP 21801	PP 11801	PP 011801	PP 021801	



# SUTURE NEEDLES

(봉합침)



**AILEE**  
SUTURES



# SUTURE NEEDLES



**AILEE**  
SUTURES



# SUTURE NEEDLES

## NEEDLES TYPE



TAPER POINT



BLUNT POINT



TAPER CUT






CUTTING EDGE



SPATULA DIAMOND



TROCAR POINT



REGULAR EYE



FRENCH EYE



FRENCH SPRING EYE



EYELESS



각막용 [Corneal]



3/8

Size No	1	2	3	4	5
Length (mm)	5	6	7	8	9


A126	3/8	Taper Point	Regular
B126	3/8	Taper Point	French
C126	3/8	Taper Point	French Spring

A125	1/2	Taper Point	Regular
B125	1/2	Taper Point	French
C125	1/2	Taper Point	French Spring

A128	3/8	Cutting	Regular
B128	3/8	Cutting	French
C128	3/8	Cutting	French Spring

A127	1/2	Cutting	Regular
B127	1/2	Cutting	French
C127	1/2	Cutting	French Spring

안과용 [Ophthalmic]



3/8

Size No	00	0	1	2	3
Length (mm)	11	13	15	18	21


A122	3/8	Taper Point	Regular
B122	3/8	Taper Point	French
C122	3/8	Taper Point	French Spring

A121	1/2	Taper Point	Regular
B121	1/2	Taper Point	French
C121	1/2	Taper Point	French Spring

A124	3/8	Cutting	Regular
B124	3/8	Cutting	French
C124	3/8	Cutting	French Spring

A123	1/2	Cutting	Regular
B123	1/2	Cutting	French
C123	1/2	Cutting	French Spring

혈관용 [Vascular]



Straight


Size No	1	2	3
Length (mm)	15	20	25

A21	Straight	Taper Point	Regular
B21	Straight	Taper Point	French
C21	Straight	Taper Point	French Spring

A23	Straight	Taper Cut	Regular
B23	Straight	Taper Cut	French
C23	Straight	Taper Cut	French Spring

# SUTURE NEEDLES

## 혈관용 [Vascular]




1/2

Size No	0	1	2	3
Length (mm)	11	14	19	24

A22	3/8	Taper Point	Regular
B22	3/8	Taper Point	French
C22	3/8	Taper Point	French Spring

A24	1/2	Taper Point	Regular
B24	1/2	Taper Point	French
C24	1/2	Taper Point	French Spring


## 치과용 [Dental]



Fishing Hook

Size No	1	2	3
Length (mm)	12	15	18

A150	Fishing Hook	Cutting	Regular
B150	Fishing Hook	Cutting	French
C150	Fishing Hook	Cutting	French Spring




1/2

Size No	1	2	3
Length (mm)	12	15	18

A151	3/8	Cutting	Regular
B151	3/8	Cutting	French
C151	3/8	Cutting	French Spring

A152	1/2	Cutting	Regular
B152	1/2	Cutting	French
C152	1/2	Cutting	French Spring

## 피부용 [Skin]




3/8

Size No	1	2	3	4	5	6	7
Length (mm)	17	21	25	29	33	37	41

A104	3/8	Cutting	Regular
B104	3/8	Cutting	French
C104	3/8	Cutting	French Spring

A103	1/2	Cutting	Regular
B103	1/2	Cutting	French
C103	1/2	Cutting	French Spring


외과용 [General Closure]



Straight

Size No	1	2	3	4	5	6	7	8	9	10
Length (mm)	30	35	40	45	50	55	60	65	70	75

A162	Straight	Cutting	Regular
B162	Straight	Cutting	French
C162	Straight	Cutting	French Spring




Half Curved

Size No	0	1	2	3	4	5	6
Length (mm)	25	30	40	50	55	60	70

A108	Half Curved	Cutting	Regular
B108	Half Curved	Cutting	French
C108	Half Curved	Cutting	French Spring

# SUTURE NEEDLES

## 외과용 [General Closure]




Size No	1	2	3	4	5	6	7	8	9	10
Length (mm)	17	21	25	29	33	37	41	45	49	53

A102	3/8	Cutting	Regular
B102	3/8	Cutting	French
C102	3/8	Cutting	French Spring

A106	1/2	Cutting	Regular
B106	1/2	Cutting	French
C106	1/2	Cutting	French Spring


## 내장용 [Gastro-Intestinal]



Size No	1	2	3	4	5	6	7	8
Length (mm)	34	42	50	58	66	74	82	90

A134	Straight	Taper Point	Regular
B134	Straight	Taper Point	French
C134	Straight	Taper Point	French Spring

내장용 [Gastro-Intestinal]




1/2

Size No	1	2	3	4	5	6	7	8
Length (mm)	18	23	28	33	38	43	48	53

A135	3/8	Taper Point	Regular
B135	3/8	Taper Point	French
C135	3/8	Taper Point	French Spring

A136	1/2	Taper Point	Regular
B136	1/2	Taper Point	French
C136	1/2	Taper Point	French Spring

내복용 [Abdominal]



1/2


Size No	1	2	3	4	5	6	7	8	9	10
Length (mm)	17	21	25	29	33	37	41	45	49	53

A113	3/8	Taper Point	Regular
B113	3/8	Taper Point	French
C113	3/8	Taper Point	French Spring

A114	1/2	Taper Point	Regular
B114	1/2	Taper Point	French
C114	1/2	Taper Point	French Spring

# SUTURE NEEDLES

## 부인과용 [Ob Gyn]




1/2

Size No	1	2	3	4	5	6
Length (mm)	31	36	41	46	51	56

A143	3/8	Taper Point	Regular
B143	3/8	Taper Point	French
C143	3/8	Taper Point	French Spring

A144	1/2	Taper Point	Regular
B144	1/2	Taper Point	French
C144	1/2	Taper Point	French Spring




1/2

Size No	0	1	2	3
Length (mm)	30	45	60	75

A148	1/2	Blunt Point	Regular
B148	1/2	Blunt Point	French
C148	1/2	Blunt Point	French Spring

부인과용 [Ob Gyn]




1/2

Size No	1	2	3	4	5	6
Length (mm)	31	36	41	46	51	56

A141	3/8	Cutting	Regular
B141	3/8	Cutting	French
C141	3/8	Cutting	French Spring

A142	1/2	Cutting	Regular
B142	1/2	Cutting	French
C142	1/2	Cutting	French Spring

근육용 [Fascia]



3/8

Size No	1	2	3	4	5	6
Length (mm)	25	32	39	46	53	60

A145	3/8	Cutting	Regular
B145	3/8	Cutting	French
C145	3/8	Cutting	French Spring


A146	1/2	Cutting	Regular
B146	1/2	Cutting	French
C146	1/2	Cutting	French Spring

A149	3/8	Taper Point	Regular
B149	3/8	Taper Point	French
C149	3/8	Taper Point	French Spring

A147	1/2	Taper Point	Regular
B147	1/2	Taper Point	French
C147	1/2	Taper Point	French Spring

# SUTURE NEEDLES


## 해부용 [Retention]



Half Curved

Size No	1	2	3	4	5
Length (mm)	70	82	100	115	130

A 118	Half Curved	Cutting	Regular
B 118	Half Curved	Cutting	French
C 118	Half Curved	Cutting	French Spring



3/8

Size No	1	2	3	4	5	6
Length (mm)	60	66	73	80	87	93

A 116	3/8	Cutting	Regular
B 116	3/8	Cutting	French
C 116	3/8	Cutting	French Spring

A 115	1/2	Cutting	Regular
B 115	1/2	Cutting	French
C 115	1/2	Cutting	French Spring



# SURGICAL BLADES





# SURGICAL BLADES



Use Handle No.3



Use Handle No.4




## Characteristics

1. Sterilized by gamma radiation
2. Made of special stainless steel
3. Grinded keenly with newest auto-grinding machine
4. The best quality guaranteed by special heat-treatment
5. Packed with auto-packing machine
6. International standard goods



# STERILE SURGICAL BLADES



**How to handle goods**

For opening aluminium-foil package, grasp both end marked 'open' and pull right and left



**Aluminium foil packed**

Number and shape of blades are imprinted on the box.  
For taking out the single peel pack, press the perforated flap at the bottom of 100 blades box.



**Packing unit**

1 packet of 100 blades.



## AILEE COMPANY LIMITED

W. [www.ailee.kr](http://www.ailee.kr)

E. [ailee1963@ailee.kr](mailto:ailee1963@ailee.kr)

---

Head  
Office

28, Gaya-daero 196beon-gil Sasang-gu, Busan, South Korea  
T. +82)51-322-6161~5 F. +82)51-323-1543

Seoul  
Office

733, Nonhyeon-ro, Gangnam-gu, Seoul, South Korea  
T. +82)2-543-6163 F. +82)51-543-6164

Chicago  
Office

475 N. Martingale Rd., Suite 710, Schaumburg IL 60173, USA

AILEE Product catalog

\* 이 Catalogue는 주식회사 아이리의 동의 없이 수정, 변경, 복사 및 배포 할 수 없습니다.